



FLUSSO™

# FLUSSO™ TFI

To Enable a Safer Circuit Disconnect



[www.keepthepeep.com](http://www.keepthepeep.com)



# FLUSSO™ TFI

Flusso™ TFI allows for temporary interruption of gas flow to maintain a patient's positive end expiratory pressure (PEEP) and reduces potential exhaled aerosols during ventilator circuit disconnect.

- Lightweight
- Clear housing permits visualization of function
- Low deadspace (1.86 cc)
- Compact
- Minimal effort to activate
- Safety tabs to prevent over compression



## ORDERING INFORMATION

Item #	Description	Qty
420-1979	Flusso™ TFI	25/box

**Maintain PEEP • Reduce Exhaled Aerosols • Lung Protective Ventilation**



McArthur Medical Sales Inc  
1846 5th Concession Rd W  
Rockton, ON L0R 1X0

1-800-996-6674  
mmsi@mcarthurmedical.com  
www.mcarthurmedical.com



[www.keepthepeep.com](http://www.keepthepeep.com)